

# Say hello to SNAP<sup>™</sup> Therapy.

The SNAP<sup>™</sup> Therapy System combines the simplicity of advanced wound dressings with the proven benefits<sup>1</sup> of negative pressure therapy in a discreet design that won't get noticed.



**SMART** FOR CLINICIANS

**SMART** FOR PATIENTS

**SMART** FOR FACILITIES

# SNAP™

Therapy System

**SILENT:** Allows patients to continue with activities of daily living<sup>1</sup>

**DISCREET:** Pocket-sized, lightweight design can be hidden under clothes

**CONVENIENT:** Disposable, single-use components are available off-the-shelf for easy access

## CLINICAL EVIDENCE

### Armstrong<sup>1</sup>


In a multicenter RCT, 132 patients with lower extremity diabetic and venous wounds were enrolled in the study. 118 patients were treated either with SNAP™ Therapy (n=59) or V.A.C.® Therapy (n=56), with 115 patients completing the study.

- Patients were treated for up to 16 weeks or complete wound closure.
- Primary end point analysis of wound size reduction found that SNAP™ Therapy treated subjects demonstrated non-inferiority to V.A.C.® Therapy subjects at 4, 8, 12 and 16 weeks ( $p=0.0030, 0.0130, 0.0051$  and  $0.0044$ , respectively).
- The study indicated that the effect of the SNAP™ Therapy System was not significantly different than that of the V.A.C.® Therapy System in promoting complete wound closure in the population studied ( $p=0.9620$ ).
- SNAP™ Therapy patients reported less interruption of activities on daily living compared to V.A.C.® Therapy patients. However, pain associated with treatment was not significantly different between treatment groups.
- Other benefits noted by the authors were shorter time to dressing application and ease of use.
- However, despite randomization, the initial wound size was significantly greater in the V.A.C.® Therapy patients than in the SNAP™ Therapy patients (mean of  $9.95\text{cm}^2$  vs  $5.37\text{cm}^2$ ;  $p=0.0093$ )

### Marston<sup>2</sup>


In a multicenter RCT, 40 patients with venous leg ulcers were treated either with SNAP™ Therapy (n=19) or V.A.C.® Therapy (n=21)

- Patients were evaluated for 16 weeks or complete wound closure.
- Primary end point analysis of wound size reduction found that SNAP™ Therapy treated subjects significantly greater wound size reduction than in V.A.C.® Therapy subjects at 4, 8, 12 and 16 weeks ( $p=0.0039, 0.0086, 0.0002$ , and  $0.0005$ , respectively).
- 53% of SNAP™ Therapy patients achieved 50% wound closure at 30 days compared to 24% of V.A.C.® Therapy patients.
- However, despite randomization, the initial wound size was significantly greater in the V.A.C.® Therapy patients than in the SNAP™ Therapy patients (mean of  $11.6\text{cm}^2$  vs  $4.49\text{cm}^2$ ).




Clean wound and apply BLUE foam.

Place Hydrocolloid Dressing over wound and seal.



Connect tubing from Hydrocolloid Dressing to SNAP™ Cartridge.

Remove Activation Key from Cartridge to initiate negative pressure therapy.



### SNAP™ Therapy Cartridge

Catalog Number	Pressure	Capacity
SNPA125US	-125mmHg	60ml

### SNAP™ Advanced Dressing Kit

Catalog Number	Size	Interface
SKTF10X10	10cm x 10cm	Foam
SKTF15X15	15cm x 15cm	Foam

### SNAP™ Bridge Dressing Kit

Catalog Number	Size	Interface
BKTF14X11	14cm x 11cm	Foam
BKTF14X11S	14cm x 11cm with SNAP™ SecurRing™ Hydrocolloid	Foam

### SNAP™ SecurRing™ Hydrocolloid

Catalog Number	Size
SRNG10	2" diameter

### SNAP™ Therapy Strap

Catalog Number	Size
STPAS	Small 18"
STPAM	Medium 21"
STPAL	Large 24"

To order product or for more information, call **800-275-4524** or visit **acelity.com**

### REFERENCES

1. Armstrong DG, Marston WA, Reyzelman AM, Kirsner RS. Comparative effectiveness of mechanically and electrically powered negative pressure wound therapy devices: a multicenter randomized controlled trial. *Wound Rep Reg.* 2012; 20(3):332-341
2. Marston WA, Armstrong DG, Reyzelman AM, Kirsner RS. A Multicenter Randomized Controlled Trial Comparing Treatment of Venous Leg Ulcers Using Mechanically Versus Electrically Powered Negative Pressure Wound Therapy. *Advances in Wound Care.* 2015;4(2):75-82. doi:10.1089/wound.2014.0575.



**NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult a physician and product instructions for use prior to application. Rx only.**

Copyright 2015, 2016 KCI Licensing, Inc. All rights reserved. All trademarks designated herein are proprietary to KCI Licensing, Inc., its affiliates and/or licensors. DSL#15-0750.US • LIT#29-A-261 (Rev. 9/16)

