

Better care for every person — all skin, every wound, every time. It's at the heart of what you do every day, and we're on a mission to help. As a leader in advanced wound care and skin integrity, we combine our science with your relentless commitment to care. Together, we can elevate the critical role skin plays in the well-being of every person, and move care forward to help patients heal from treatment and get back to their lives.

Protect and Manage Skin

3M™ Cavilon™ Advanced Skin Protectant

3M™ Cavilon™ Advanced Skin Protectant utilizes revolutionary polymer-cyanoacrylate technology to create a highly durable, ultra-thin barrier that protects damaged or intact skin from caustic, corrosive body fluids such as liquid stool and gastric fluid.



Item No.	Size	Each/Case	HCPCS Code
5050	2.7 mL applicator	20	A6250
5051	0.7 mL applicator	20	—

3M™ Cavilon™ No Sting Barrier Film

3M™ Cavilon™ No Sting Barrier Film is the original alcohol-free barrier film. It is a gentle, effective and CHG-compatible solution for routine skin protection and is supported by more than 80 pieces of evidence.



Item No.	Size	Each/Case	HCPCS Code
3343	1 mL wand	25 items/box, 4 boxes/case	Skin Care: A6250 Ostomy: A5120
3344	1 mL wipe	30 items/box, 4 boxes/case	Skin Care: A6250 Ostomy: A5120
3345	3 mL wand	25 items/box, 25 boxes/case	Skin Care: A6250 Ostomy: A5120
3346	28 mL spray bottle	12 items/box, 1 box/case	Skin Care: A6250 Ostomy: A5120

3M™ Cavilon™ No-Rinse Skin Cleanser

3M™ Cavilon™ No-Rinse Skin Cleanser is a no-rinse cleanser that can be used for incontinence and general skin cleansing needs.



Item No.	Size	Each/Case	HCPCS Code
3380	8 oz. bottle	12	A6250

Protect and Manage Skin (cont.)

3M™ Cavilon™ Durable Barrier Cream

Formulated with 3M polymer technology, 3M™ Cavilon™ Durable Barrier Cream is a concentrated, fragrance-free moisturizing skin barrier that resists wash-off and allows adhesives to adhere. It's ideal for moisturization of at-risk and/or severely dry skin. It also provides a barrier which protects skin at low risk for breaking down.



Item No.	Size	Each/Case	HCPCS Code
3353	2 g sachet	20 items/box, 12 boxes/case	A6250
3354	1 oz. tube	48 items/box, 1 box/case	A6250
3355	3.25 oz. tube	12 items/box, 1 box/case	A6250

3M™ Cavilon™ Antifungal Cream

3M™ Cavilon™ Antifungal Cream effectively provides relief from redness, irritation, scaling, itching, discomfort and burning caused by fungal skin infections.



Item No.	Size	Each/Case	HCPCS Code
3389	2 oz. tube	24	A6250
3390	5 oz. tube	12	A6250

3M™ Cavilon™ Extra Dry Skin Cream

This versatile skin cream can be used anywhere on the body to relieve extremely dry skin.



Item No.	Size	Each/Case	HCPCS Code
3386	4 oz. tube	12	A6250

3M™ Cavilon™ Moisturizing Hand Lotion

3M™ Cavilon™ Moisturizing Hand Lotion is a skin friendly, everyday moisturizer that's hypoallergenic and fragrance-free.



Item No.	Size	Each/Case	HCPCS Code
9205	16 oz. pump bottle	12	A6250
9215	2 oz. personal size bottle	48	A6250

Manage Biofilm/Bioburden

TOUCHLESS CARE™ Spray Family

TOUCHLESS CARE™ Sprays are designed with formulas that protect skin, maintain moisture balance and allow for a convenient spray application.



Item No.	Size	Description	Each/Case
62402	2.0 oz. bottle	Zinc Oxide Protectant Spray	24
62404	4.5 oz. bottle	Zinc Oxide Protectant Spray	24
82402	2.0 oz. bottle	Antifungal Spray	24

KERRACEL™ Ag Gelling Fiber Dressing

KERRACEL™ Ag Dressing is a gelling fiber silver dressing that uses the combination of an enhanced carboxymethyl cellulose (CMC) dressing and the advanced silver technology — Ag OXYSALTS™ Technology* — to create the optimal wound healing environment.




Item No.	Size	Each/Case	HCPCS Code
CWL1157	2" x 2"	10	A6196
CWL1158	4" x 5"	10	A6197
CWL1159	6" x 6"	5	A6197
CWL1160	8" x 12"	5	A6198
CWL1168	0.75" x 12" ribbon	5	A6199
CWL1162	0.75" x 18" ribbon	5	A6199
CWL1163	16" x 16"	5	—

*Ag OXYSALTS™ Technology is a trademark of Argent Limited Partnership.

SILVERCEL™ NON-ADHERENT Antimicrobial Alginate Dressing with EASYLIFT™ Precision Film Technology

SILVERCEL™ NON-ADHERENT Dressing is a non-woven pad composed of alginate, carboxymethylcellulose (CMC) and silver coated nylon fibers, with a non-adherent wound contact layer.



Item No.	Size	Each/Case	HCPCS Code
900112	1" x 12"	5 items/box, 5 boxes/case	A6199
900202	2" x 2"	10 items/box, 5 boxes/case	A6196
900404	4-1/4" x 4-1/4"	10 items/box, 5 boxes/case	A6197
900408	4" x 8"	5 items/box, 5 boxes/case	A6197

KERRACONTACT™ Ag Dressing

KERRACONTACT™ Ag Dressing is a silver dressing for chronic and acute wounds that utilizes Ag OXYSALTS™ Technology.* Ag OXYSALTS™ Technology is a unique type of silver that is up to 6 times more powerful than Ag¹⁺ while at the same time being gentle on patients and your budget.**




Item No.	Size	Each/Carton	HCPCS Code
PN-09-0102	2" x 2"	10	A6206
PN-09-0103	4" x 5"	10	A6207
PN-09-0104	4" x 8"	10	A6207
PN-09-0105	6" x 6"	10	A6207
PN-09-0106	6" x 9"	10	A6208
PN-09-0108	8" x 16"	10	A6208
PN-09-0116	16" x 16"	5	A6208
PN-09-0120	4" x 48"	5	A6208
PN-09-0132	16" x 32"	5	A6208

*Ag OXYSALTS™ Technology is a trademark of Argent Limited Partnership.
 **Spina C. Silver I, II, III: Chemical Characteristics, Properties, and Antimicrobial Activity. Poster presented at SAWC Spring; April 23–27, 2014; Orlando, FL.

Manage Biofilm/Bioburden (cont.)

ACTISORB™ Silver 220 Antimicrobial Binding Dressing

ACTISORB™ Silver 220 Antimicrobial Binding Dressing is an activated charcoal dressing with silver, enclosed in a non-adherent nylon sleeve.



Item No.	Size	Each/Box	HCPCS Code
650220	2.6" x 3.8"	10 items/carton, 5 cartons/box	A4649
105220	4.2" x 4.2"	10 items/carton, 5 cartons/box	A4649
190220	4.2" x 7.6"	10 items/carton, 5 cartons/box	A4649

Provide Collagen

PROMOGRAN™ Matrix Wound Dressing/PROMOGRAN PRISMA™ Matrix

PROMOGRAN™ Matrix is comprised of 45% oxidized regenerated cellulose (ORC) and 55% collagen in a sterile, freeze dried composite.

PROMOGRAN PRISMA™ Matrix is comprised of 44% oxidized regenerated cellulose (ORC), 55% collagen and 1% silver-ORC in a sterile, freeze dried composite.


PROMOGRAN PRISMA™ Matrix contains a compound of silver and ORC and has nominally twice the amount of Collagen/ORC material as the same size of PROMOGRAN™ Matrix.



Item No.	Size	Each/Box	HCPCS Code
PROMOGRAN™ Matrix			
PG004	4.34 sq. in. hexagon	10 items/carton, 4 cartons/box	A6021
PG019	19.1 sq. in. hexagon	10 items/carton, 4 cartons/box	A6022
PROMOGRAN PRISMA™ Matrix			
MA028	4.34 sq. in. hexagon	10 items/carton, 4 cartons/box	A6021
MA123	19.1 sq. in. hexagon	10 items/carton, 4 cartons/box	A6022
MA032	3/8" x 3/8" x 12-5/8" rope	6 items/carton, 4 cartons/box	A6024

FIBRACOL™ Plus Collagen Wound Dressing with Alginate

FIBRACOL™ Plus Dressing with Alginate is a soft, absorbent, and conformable wound dressing composed of 90% collagen and 10% calcium alginate.




Item No.	Size	Each/Box	HCPCS Code
2981	2" x 2"	12 items/carton, 6 cartons/box	A6021
2982	4" x 4-3/8"	12 items/carton, 6 cartons/box	A6022
2983	4" x 8-3/4"	6 items/carton, 6 cartons/box	A6022
2984	3/8" x 3/8" x 15-3/4"	6 items/carton, 6 cartons/box	A6024

Minimize Adherence

ADAPTIC TOUCH™ Non-Adhering Silicone Dressing

ADAPTIC TOUCH™ Dressing is a flexible, open-mesh silicone wound contact layer designed to optimize fluid management and minimize pain at dressing change.



Item No.	Size	Each/Box	HCPCS Code
500501	2" x 3"	10 items/carton, 5 cartons/box	A6206
500502	3" x 4-1/4"	10 items/carton, 5 cartons/box	A6206
500503	5" x 6"	10 items/carton, 5 cartons/box	A6207
500504	8" x 12-3/4"	5 cartons/box	A6208

ADAPTIC™ Digit Non-Adhering Dressing


ADAPTIC™ Digit Dressing is made of a unique, tacky, silicone-coated fabric covered with a tubular bandage. It is specifically designed for dressing fingers and toes.



Item No.	Size	Each/Box
ADAPTIC™ Digit Dressing		
MAD003U	Small 3/4"	10
MAD013U	Medium 1"	10
MAD023U	Large 1-1/8"	10
MAD042U	XL 1-1/4"	10
ADAPTIC™ Digit Toe		
MAD062U	Large 1-1/8"	10

ADAPTIC™ Non-Adhering Dressing

ADAPTIC™ Dressing is a primary dressing made of knitted cellulose acetate mesh, impregnated with a specially formulated petrolatum emulsion.



Item No.	Size	Each/Box	HCPCS Code
2012	3" x 3"	50 items/carton, 12 cartons/box	A6222
2013	3" x 8"	36 x 3 items/carton, 6 cartons/box	A6223
2014	3" x 16"	36 items/carton, 6 cartons/box	A6223
2015	3" x 8"	24 items/carton, 6 cartons/box	A6223
2018	3" x 60" roll	10 rolls/box	A6266
2019	5" x 9" sheet	12 sheets/carton, 6 cartons/box	A6223

Manage Exudate

3M™ Tegaderm™ Transparent Film Roll

3M™ Tegaderm™ Transparent Film Roll provides maximum versatility for waterproof protection of skin and cover dressings.* It protects at-risk, intact skin to help prevent skin breakdown caused by friction and bodily fluids.



*Product is non-sterile; not recommended as a primary wound or IV dressing.

Item No.	Size	Each/Case	HCPCS Code
16002	2" x 11 yd. (5 cm x 10 m)	1 roll/box, 4 boxes/case	A6257-A6259
16004	4" x 11 yd. (10 cm x 10 m)	1 roll/box, 4 boxes/case	A6257-A6259
16006	6" x 11 yd. (15 cm x 10 m)	1 roll/box, 4 boxes/case	A6257-A6259

KERRALITE COOL™ Hydrogel Moisture Balancing Dressing

KERRALITE COOL™ Dressings are soothing, debriding and moisturizing dressings that provide the ideal environment for managing wounds with light exudate.



Item No.	Size	Each/Box	HCPCS Code
KERRALITE COOL™ Non Border Dressing			
CWL1004	2.4" x 2.4"	5	A6234
CWL1005	5" x 3.5"	5	A6235
CWL1006	7" x 5"	5	A6235
KERRALITE COOL™ Border Dressing			
CWL1007	3" x 3"	5	A6237
CWL1008	4" x 4"	5	A6237
CWL1009	6" x 6"	5	A6238

NU-GEL™ Hydrogel with Alginate

Put hydration where your patient needs it with NU-GEL™ Hydrogel. It creates a moist wound healing environment, which assists with natural autolytic debridement while the alginate component enhances its absorptive capabilities. The gel can be used to soften and hydrate eschar by facilitating rehydration of the wound.



Item No.	Size	Each/Case	HCPCS Code
MNG415U	15 g (0.52 oz.)	10 ampules/box, 4 boxes/case	A6248
MNG425U	25 g (0.88 oz.)	6 ampules/box, 4 boxes/case	A6248

3M™ Tegaderm™ Hydrogel Wound Filler

3M™ Tegaderm™ Hydrogel Wound Filler is an easy-to-apply, preservative-free, amorphous hydrogel formulated to help provide a moist wound healing environment that has been shown to enhance wound healing.



Item No.	Size	Each/Case	HCPCS Code
91110	15 g (0.52 oz.)	10 tubes/box, 10 boxes/case	A6248
91111	25 g (0.88 oz.)	10 tubes/box, 10 boxes/case	A6248

Manage Exudate (cont.)

3M™ Tegaderm™ Absorbent Clear Acrylic Dressing

3M™ Tegaderm™ Absorbent Clear Acrylic Dressing is a unique, absorbent dressing that offers a clear view of the wound or surgical incision. Offers a waterproof, transparent barrier — allowing active monitoring and early interventions without disturbing the wound or incision.



*Dimensions of entire dressing, including hydrocolloid and film backing.

Item No.	Dressing Pad Size	Overall Dressing Size	Each/Case	HCPCS Code
90800 (Oval)	1-1/2" x 2-1/4" (3.8 cm x 5.7 cm)	3" x 3-3/4" (7.6 cm x 9.5 cm)	5 dressings/box, 6 boxes/case	A6203
90801 (Oval)	2-3/8" x 3" (6.0 cm x 7.6 cm)	4-3/8" x 5" (11.1 cm x 12.7 cm)	5 dressings/box, 6 boxes/case	A6203
90802 (Square)	3-7/8" x 4" (9.8 cm x 10.1 cm)	5-7/8" x 6" (14.9 cm x 15.2 cm)	5 dressings/box, 6 boxes/case	A6203
90803 (Oval)	3-3/8" x 4-1/4" (8.5 cm x 10.7 cm)	5-5/8" x 6-1/4" (14.2 cm x 15.8 cm)	5 dressings/box, 6 boxes/case	A6203
90805 (Square)	5-7/8" x 6" (14.9 cm x 15.2 cm)	7-7/8" x 8" (20.0 cm x 20.3 cm)	5 dressings/box, 4 boxes/case	A6204
90815 (Rectangle)	1-3/4" x 11-3/4" (4.5 cm x 20.0 cm)	3-1/2" x 10" (9.0 cm x 25.0 cm)	30 dressings/box, 4 boxes/case	A6203
90817 (Rectangle)	1-3/4" x 11-3/4" (4.5 cm x 30.0 cm)	3-1/2" x 13-3/4" (9.0 cm x 35.0 cm)	30 dressings/box, 4 boxes/case	A6204

3M™ Tegaderm™ Hydrocolloid Dressing

3M™ Tegaderm™ Hydrocolloid Dressing is designed for low to moderate draining wounds. It features up to seven-day wear time.



*Dimensions of entire dressing, including hydrocolloid and film backing.

Cat. No.	Gel Pad Size	Dressing Size	Each/Case	HCPCS Code
90001 (Ovals)	2-3/4" x 3-1/2" (6.9 cm x 8.8 cm)	4" x 4-3/4"* (10.1 cm x 12 cm)	5 dressings/box, 20 boxes/case	A6237
90003 (Ovals)	4" x 4-3/4" (10.1 cm x 12.0 cm)	5-1/4" x 6"* (13.3 cm x 15.2 cm)	5 dressings/box, 12 boxes/case	A6238
90004 (Ovals)	5-1/2" x 6-3/4" (13.9 cm x 17.1 cm)	6-3/4" x 8"* (17.1 cm x 20.3 cm)	3 dressings/box, 20 boxes/case	A6238
90002 (Squares)	4" x 4" (10.1 cm x 10.1 cm)	4" x 4" (10.1 cm x 10.1 cm)	5 dressings/box, 20 boxes/case	A6234
90005 (Squares)	6" x 6" (15.2 cm x 15.2 cm)	6" x 6" (15.2 cm x 15.2 cm)	3 dressings/box, 20 boxes/case	A6235

3M™ Tegaderm™ Hydrocolloid Dressing Sacral Design

3M™ Tegaderm™ Hydrocolloid Dressing Sacral Design is shaped specifically for one-handed application to the sacral-coccyx area to allow long wear time. It's designed for low to moderate draining wounds. It features up to seven-day wear time.



*Dimensions of entire dressing, including hydrocolloid and film backing.

Item No.	Gel Pad Size	Dressing Size	Each/Case	HCPCS Code
90007	4-7/8" x 5-1/2" (12.3 cm x 15.9 cm)	6-3/4" x 6-3/8" (17.1 cm x 16.1 cm)*	6 dressings/box, 4 boxes/case	A6238

3M™ Tegaderm™ Hydrocolloid Thin Dressing

3M™ Tegaderm™ Hydrocolloid Thin Dressing is a wound dressing that consists of a hydrocolloid adhesive with an outer clear cover film. This transparent dressing is used for mild draining wounds on at-risk skin or on intact skin to protect from friction.




*Dimensions of entire dressing, including hydrocolloid and film backing.

Item No.	Gel Pad Size	Dressing Size	Each/Case	HCPCS Code
90021 (Ovals)	2-3/4" x 3-1/2" (6.9 cm x 8.8 cm)	4" x 4-3/4"* (10.1 cm x 12.0 cm)	10 dressings/box, 10 boxes/case	A6237
90023 (Ovals)	4" x 4-3/4" (10.1 cm x 12 cm)	5-1/4" x 6"* (13.3 cm x 15.2 cm)	10 dressings/box, 6 boxes/case	A6237
90022 (Squares)	4" x 4" (10.1 cm x 10.1 cm)	4" x 4" (10.1 cm x 10.1 cm)	5 dressings/box, 20 boxes/case	A6234

Manage Exudate (cont.)

3M™ Tegaderm™ Silicone Foam Dressing

3M™ Tegaderm™ Silicone Foam Dressings are an excellent choice for wound management and may help prevent skin damage as part of a comprehensive pressure ulcer prevention program.



Featuring 3M's innovative adhesive technology, they offer a significantly longer wear time than the leading competitive silicone foam dressings* while being gentle to the skin. They feature a highly conformable low-profile edge to help minimize rolling and lifting, especially in high shear locations. Unique multi-layer design that absorbs and evaporates moisture with a patented spoke delivery system enabling easy application. The dressing may be lifted and repositioned to allow for skin assessment.

*4x4 and 6x6 dressings, based on *In vivo* studies EM-13977 and EM-13978. 3M data on file.

Item No.	Product	Overall Dressing Size	Each/Case	HCPCS Code
90631	Non-bordered Dressing	4" x 4-1/4" (10 cm x 11 cm)	10 dressings/box, 4 boxes/case	A6210
90632	Non-bordered Dressing	6" x 6" (15 cm x 15 cm)	10 dressings/box, 4 boxes/case	A6210
90643	Bordered Dressing	2" x 2" (5 cm x 5 cm)	10 dressings/box, 6 boxes/case	A6212
90640	Bordered Dressing	3" x 3" (8 cm x 8 cm)	10 dressings/box, 6 boxes/case	A6212
90641	Bordered Dressing	4" x 4" (10 cm x 10 cm)	10 dressings/box, 6 boxes/case	A6212
90642	Bordered Dressing	6" x 6" (15 cm x 15 cm)	10 dressings/box, 4 boxes/case	A6213
90646	Heel and Contour	6" x 6" (15 cm x 15 cm)	5 dressings/box, 4 boxes/case	A6212
90647	Small Sacral	6" x 6-3/4" (15 cm x 17 cm)	10 dressings/box, 4 boxes/case	A6213
90648	Large Sacral	7-1/4" x 8-3/4" (18.5 cm x 22 cm)	5 dressings/box, 4 boxes/case	A6213

3M™ Tegaderm™ High Performance Foam Adhesive Dressing

3M™ Tegaderm™ High performance Foam Adhesive Dressings are designed to meet the challenges of low- to high-exuding wounds. Unique multi-layer design provides high absorbency with high breathability to reduce the risk of skin maceration. Film backing prevents strike-through and protects against external bacteria and viruses.*



*In vitro testing shows that the transparent film provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.

**The mini wrap dressing is constructed of a conformable, absorbent, polyurethane foam pad with a highly breathable, non-waterproof, film backing reinforced with soft cloth tape.

Item No.	Foam Pad Size	Overall Dressing Size	Each/Case	HCPCS Code
90610 (Square)	2" x 2" (5 cm x 5 cm)	3-1/2" x 3-1/2" (8.8 cm x 8.8 cm)	10 dressings/box, 4 boxes/case	A6212
90611 (Oval)	2-1/2" x 3" (6 cm x 7.6 cm)	4" x 4-1/2" (10 cm x 11 cm)	10 dressings/box, 4 boxes/case	A6212
90612 (Square)	4" x 4" (10 cm x 10 cm)	5-5/8" x 5-5/8" (14.3 cm x 14.3 cm)	10 dressings/box, 4 boxes/case	A6212
90613 (Oval)	4" x 4-1/2" (10 cm x 11 cm)	5-5/8" x 6-1/8" (14.3 cm x 15.6 cm)	5 dressings/box, 6 boxes/case	A6213
90614 (Mini Oval)	1-1/4" x 1-1/2" (3.1 cm x 3.8 cm)	2-3/4" x 3" (6.9 cm x 7.6 cm)	10 dressings/box, 4 boxes/case	A6212
90616 (Oval)	5-1/2" x 6-3/4" (14 cm x 17.1 cm)	7-1/2" x 8-3/4" (19 cm x 22.2 cm)	5 dressings/box, 3 boxes/case	A6213
90619 (Heel/Elbow)	3" x 3" (7.6 cm x 7.6 cm)	5-1/2" x 5-1/2" (13.9 cm x 13.9 cm)	5 dressings/box, 4 boxes/case	A6212
90615** (Mini Wrap)	1" x 1" (2.5 cm x 2.5 cm)	2-3/4" x 2-3/4" (7 cm x 7 cm)	10 dressings/box, 4 boxes/case	A6212

3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing

3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressings are designed to meet the challenges of low- to high-exuding wounds. Innovative layer technology absorbs and evaporates moisture to maintain an optimal healing environment.



Item No.	Size	Each/Case	HCPCS Code
90600	2" x 2" (5 cm x 5 cm)	10 dressings/box, 4 boxes/case	A6209
90601	4" x 4" (10 cm x 10 cm)	10 dressings/box, 4 boxes/case	A6209
90602	4" x 8" (10 cm x 20 cm)	5 dressings/box, 6 boxes/case	A6210
90603	8" x 8" (20 cm x 20 cm)	5 dressings/box, 6 boxes/case	A6211
90604	3-1/2" x 3-1/2" (8.8 cm x 8.8 cm), Fenestrated	10 dressings/box, 4 boxes/case	A6209
90605	4" x 24" roll (10 cm x 60 cm)	1 roll/package, 6 packages/case	A6211

Manage Exudate (cont.)

3M™ Tegaderm™ High Integrity and High Gelling Alginate Dressing

3M™ Tegaderm™ High Integrity and High Gelling Alginate Dressings are optimal for maintaining a moist wound healing environment. The absorbent dressing may be used on moderate to heavily exuding wounds, absorbing up to 20 times their own weight in exudate.



Item No.	Size	Each/Case	HCPCS Code
3M™ Tegaderm™ High Integrity Alginate Dressing			
90112	4" x 4" (10 cm x 10 cm)	10 dressings/box, 5 boxes case	A6196
90120	12" rope (30.4 cm)	5 dressings/box, 4 boxes case	A6199
3M™ Tegaderm™ High Gelling Alginate Dressing			
90212	4" x 4" (10 cm x 10 cm)	10 dressings/box, 5 boxes case	A6196
90220	12" rope (30.4 cm)	5 dressings/box, 4 boxes case	A6199

KERRACEL™ Gelling Fiber Dressing


KERRACEL™ Dressing is a versatile carboxymethyl cellulose (CMC) primary dressing for moderate exuding wounds.




Item No.	Size	Each/Box	HCPCS Code
CWL1032	2" x 2"	10	A6196
CWL1033	4" x 5"	10	A6197
CWL1034	6" x 6"	5	A6197
CWL1035	1" x 18" ribbon	5	A6199

KERRAMAX CARE™ Super-Absorbent Dressing

KERRAMAX CARE™ Super-Absorbent Dressing is a primary or secondary dressing, suitable for chronic and acute wounds and helps to improve wound healing.



KERRAMAX CARE™
Super-Absorbent Dressing



KERRAMAX CARE™
Gentle Border Dressing

Item No.	Size	Each/Box	HCPCS Code
KERRAMAX CARE™ Super-Absorbent Dressing			
PRD500-025	2" x 2"	10	A6196
PRD500-050	4" x 4"	10	A6196
PRD500-100	5" x 6"	10	A6197
PRD500-120	4" x 9"	10	A6197
PRD500-240	8" x 9"	10	A6197
PRD500-380	8" x 12"	5	A6198
PRD500-600	8" x 20"	5	A6198
PRD500-300	8" x 9" multisite	5	A6197
KERRAMAX CARE™ Gentle Border Dressing			
PRD500-1174	4" x 4"	10	A6196
PRD500-1175	6" x 6"	10	A6196
PRD500-1176	8" x 8"	10	A6197
PRD500-1177	6" x 10"	10	A6197

Provide Compression

3M™ Coban™ 2 and 3M™ Coban™ 2 Lite Two-Layer Compression Systems

3M™ Coban™ 2 Two-Layer Compression Systems are easy to use and provide high compression (35–40 mmHg) for patients with ABPI greater than or equal to 0.8, and reduced compression (25–30 mmHg) for patients with ABPI greater than or equal to 0.5. Each system is supplied as a kit that includes two rolls: a Comfort Layer roll and a Compression Layer roll. The thin, lightweight and breathable sleeve allows patients to wear their own footwear. 3M™ Coban™ 2 Compression Systems are comfortable to wear and provide sustained therapeutic compression for up to seven days.



Application/Use	Item No.	Comfort Layer Size		Each/Case	HCPCS Code
		Compression Layer Size			
Below the knee	2094N	4" x 2.9 yd. (10 cm x 2.7 m)		8	A6441
		4" x 5.1 yd. (10 cm x 4.7 m) stretched			A6452
Below the knee (ABPI > 0.5)	2794N	4" x 2.9 yd. (10 cm x 2.7 m)		8	A6441
		4" x 5.1 yd. (10 cm x 4.7 m) stretched			A452
Below the knee (Lg. circumference)	2094XL	4" x 3.8 yd. (10 cm x 3.5 m)		8	A6441
		4" x 6.3 yd. (10 cm x 5.8 m) stretched			A6452
Above the knee	20096	6" x 3.8 yd. (15 cm x 3.5 m)		8	—
		6" x 4.9 yd. (15 cm x 4.5 m) stretched			—
Toe boot	2092	2" x 1.3 yd. (5 cm x 1.2 m)		8	—
		2" x 3 yd. (5 cm x 2.7 m) stretched			—

Miscellaneous Products

3M™ Tegaderm™ +Pad Dressing with Non-Adherent Pad

3M™ Tegaderm™ +Pad Dressing with Non-Adherent Pad is an all-in-one dressing that provides a sterile, waterproof, viral* and bacterial barrier. These dressings consist of a non-adherent absorbent pad bonded to a thin film transparent dressing.



**In vitro* testing shows that the transparent film provides a viral barrier from viruses 27nm in diameter or larger while the dressing remains intact without leakage.

Item No.	Dressing Pad Size	Overall Dressing Size	Each/Case	HCPCS Code
3582	1" x 1-1/2" (2.5 cm x 4 cm)	2" x 2-3/4" (5 cm x 7 cm)	50 dressings/box, 4 boxes/case	A6203
3584	1" x 2-3/8" (2.5 cm x 6 cm)	2-3/8" x 4" (6 cm x 10 cm)	50 dressings/box, 4 boxes/case	A6203
3586	1-3/4" x 2-3/8" (4.5 cm x 6 cm)	3-1/2" x 4" (9 cm x 10 cm)	25 dressings/box, 4 boxes/case	A6203
3587 (Oval)	1-3/4" x 2-3/8" (4.5 cm x 6 cm)	3-1/2" x 4-1/8" (9 cm x 10.5 cm)	25 dressings/box, 4 boxes/case	A6203
3588	4" x 4" (10 cm x 10 cm)	6" x 6" (15 cm x 15 cm)	25 dressings/box, 4 boxes/case	A6203
3589	1-3/4" x 6" (4.5 cm x 15 cm)	3-1/2" x 6" (9 cm x 15 cm)	25 dressings/box, 4 boxes/case	A6203
3590	1-3/4" x 6" (4.5 cm x 15 cm)	3-1/2" x 8" (9 cm x 20 cm)	25 dressings/box, 4 boxes/case	A6203
3591	1-3/4" x 8" (4.5 cm x 20 cm)	3-1/2" x 10" (9 cm x 25 cm)	25 dressings/box, 4 boxes/case	A6203
3593	1-3/4" x 11-3/4" (4.5 cm x 30 cm)	3-1/2" 13-3/4" (9 cm x 35 cm)	25 dressings/box, 4 boxes/case	A6204

Miscellaneous Products (cont.)

3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing

3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressings are gentle, comfortable and absorbent. These dressings are porous and breathable for increased patient comfort — allowing air in and moisture out. They are easy to apply and gentle to the skin.



Item No.	Dressing Pad Size	Overall Dressing Size	Each/Case	HCPCS Code
3562	1" x 1-1/2" (2.5 cm x 4 cm)	2" x 2-3/4" (5 cm x 7 cm)	50 dressings/box, 4 boxes/case	A6219
3564	1" x 2-3/8" (2.5 cm x 6 cm)	2-3/8" x 4" (6 cm x 10 cm)	50 dressings/box, 4 boxes/case	A6219
3566	1-3/4" x 2-3/8" (4.5 cm x 6 cm)	3-1/2" x 4" (9 cm x 10 cm)	25 dressings/box, 4 boxes/case	A6219
3568	4" x 4" (10 cm x 10 cm)	6" x 6" (15 cm x 15 cm)	25 dressings/box, 4 boxes/case	A6219
3569	1-3/4" x 4" (4.5 cm x 10 cm)	3-1/2" x 6" (9 cm x 15 cm)	25 dressings/box, 4 boxes/case	A6219
3570	1-3/4" x 6" (4.5 cm x 15 cm)	3-1/2" x 8" (9 cm x 20 cm)	25 dressings/box, 4 boxes/case	A6219
3571	1-3/4" x 8" (4.5 cm x 20 cm)	3-1/2" x 10" (9 cm x 25 cm)	25 dressings/box, 4 boxes/case	A6219
3573	1-3/4" x 11-3/4" (4.5 cm x 30 cm)	3-1/2" 13-3/4" (9 cm x 35 cm)	25 dressings/box, 4 boxes/case	A6220

3M™ Steri-Strip™ Antimicrobial Skin Closure

3M™ Steri-Strip™ Antimicrobial Skin Closures are small but effective way for you to help reduce the risk of SSIs on your patients and help protect your bottom line. The iodophor adhesive has broad spectrum activity against the most common pathogens causing SSIs.



Item No.	Size	Each/Case
A1841	1/4" x 3" (6 mm x 75 mm)	3 strips/envelope, 50 envelopes/box, 4 boxes/case
A1847	1/2" x 4" (12 mm x 100 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case
A1848	1" x 5" (25 mm x 125 mm)	4 strips/envelope, 25 envelopes/box, 4 boxes/case

3M™ Steri-Strip™ Reinforced Skin Closure

3M™ Steri-Strip™ Reinforced Skin Closures are the standard for wound closure and wound support. They're sterile, comfortable, breathable adhesive strips reinforced with polymer filaments, for strong, reliable closure of incisions and skin lacerations.



Item No.	Size	Each/Case
R1540	1/8" x 3" (3 mm x 75 mm)	5 strips/envelope, 50 envelopes/box, 4 boxes/case
R1541	1/4" x 3" (6 mm x 75 mm)	3 strips/envelope, 50 envelopes/box, 4 boxes/case
R1542	1/4" x 1-1/2" (6 mm x 38 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case
R1546	1/4" x 4" (6 mm x 100 mm)	10 strips/envelope, 50 envelopes/box, 4 boxes/case
R1547	1/2" x 4" (12 mm x 100 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case
R1548	1" x 5" (25 mm x 125 mm)	4 strips/envelope, 25 envelopes/box, 4 boxes/case
R1549	1/2" x 2" (12 mm x 50 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case

3M™ Steri-Strip™ Wound Closure System

The 3M™ Steri-Strip™ Wound Closure System combines the suture and staple-free wound closure of our Steri-Strip™ products with the bacterial and viral* protection of 3M™ Tegaderm™ Transparent film dressings.



Item No.	Dressing Size	Steri-Strip Size	Each/Case
W8512	2-3/8" x 1-7/8" (60 mm x 47 mm)	1-7/8" x 1/2" (47 mm x 12 mm)	1 each/envelope, 25 envelopes/box, 4 boxes/case
W8514	2-3/8" x 1-7/8" (60 mm x 47 mm)	1-7/8" x 1/2" (47 mm x 12 mm)	3 each/envelope, 25 envelopes/box, 4 boxes/case
W8516	2-3/8" x 1-7/8" (60 mm x 47 mm)	1-7/8" x 1/2" (47 mm x 12 mm)	5 each/envelope, 25 envelopes/box, 4 boxes/case

*In vitro testing shows that the transparent film provides a viral barrier from viruses 27nm in diameter or larger while the dressing remains intact without leakage.

Miscellaneous Products (cont.)

3M™ Steri-Strip™ Blend Tone Skin Closure

Deliver wound support with a less noticeable look using 3M™ Steri-Strip™ Blend Tone Skin Closures. Their light color blends with many skin tones for a cosmetically pleasing appearance. They're sterile, breathable and comfortable to wear.



Item No.	Size	Each/Case
B1551	1/4" x 3" (6 mm x 75 mm)	3 strips/envelope, 50 envelopes/box, 4 boxes/case
B1557	1/2" x 4" (12 mm x 100 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case
B1559	1/2" x 2" (12 mm x 50 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case

3M™ Steri-Strip™ Elastic Skin Closure

Designed to stretch and accommodate post-surgical edema, body movement and contours, 3M™ Steri-Strip™ Elastic Skin Closures are a sterile and breathable non-invasive skin closure solution. Their soft, non-woven backing make them comfortable for patients to wear.



Item No.	Size	Each/Case
E4541	1/4" x 3" (6 mm x 75 mm)	3 strips/envelope, 50 envelopes/box, 4 boxes/case
E4542	1/4" x 1-1/2" (6 mm x 38 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case
E4547	1/2" x 4" (12 mm x 100 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case
E4548	1" x 5" (25 mm x 125 mm)	4 strips/envelope, 25 envelopes/box, 4 boxes/case
E4549	1/2" x 2" (12 mm x 50 mm)	6 strips/envelope, 50 envelopes/box, 4 boxes/case

3M™ Steri-Strip™ Compound Benzoin Tincture

The 3M™ Steri-Strip™ Compound Benzoin Tincture is an alcohol solution of compound benzoin that is intended for use as an adhesive adjunct for the enhancement of surgical tapes or wound closure strips.* The sterile tray packaging allows aseptic delivery.



Item No.	Description	Each/Case
C1544	2/3 cc	40 items/box, 4 boxes/case

*3M data on file: EM-05-136253.

3M™ Wound Cleanser

3M™ Wound Cleanser is a non-antimicrobial, no-rinse liquid with a zinc-nutrient formulation that provides a moist, mildly acidic environment conducive to wound healing.



Item No.	Description	Each/Case	HCPCS Code
91101	4 oz. squeeze bottle (114 mL)	12	A6250



3M Company
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

1-800-228-3957
3M.com/Medical

For clinical questions, please contact your local 3M Account Manager. To connect with us on reimbursement questions, please email us at Reimbursementeducation@mmm.com or call 1-800-668-6812.

NOTE: Specific instructions, contraindications, warnings, precautions, and safety information exist for these products and therapies, some of which may be Rx only. Please consult a clinician and product instructions for use prior to application.

© 2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.
US_70-2011-8092-7 PRA-PM-US-02861 (02/21)