


The SNAP™ Therapy System Makes a Difference in Patients' Quality of Life, Without Being Visible.

The SNAP™ Therapy System is indicated for patients who would benefit from wound management via the application of negative pressure, particularly as the device may promote wound healing through the removal of small amounts of excess exudate, infectious material, and tissue debris.

The SNAP™ Therapy System is indicated for removal of exudate from:



SMART FOR CLINICIANS

SMART FOR PATIENTS

SMART FOR FACILITIES

For more information, call **800-275-4524** or visit **acelity.com**