

KCI is proud to offer an unparalleled portfolio of advanced wound dressings and skin care products to help manage your patient's wound care needs. With a legacy of continuous innovation in wound care, KCI has a range of products that may help provide antimicrobial benefit, exudate management, compression, skin protection, minimize adherence, or anything in between- regardless of what your patient needs may be - we have you covered!

SKIN CARE



TOUCHLESS CARE™ Spray Barriers

TOUCHLESS CARE™ Sprays are designed with advanced concentrated formulas that optimize the therapeutic effect¹, and help maintain a moisture balance.² These formulations allow for the convenience of a spray application.

TOUCHLESS CARE™ Spray Barriers

Item code	Size	Description	Bottles/Box
62402	2.0oz	Zinc Oxide Protectant Spray	24
62404	4.5oz	Zinc Oxide Protectant Spray	24
72402	2.0oz	Clear Protectant Spray	24
82402	2.0oz	Antifungal Spray	24

ANTIMICROBIAL DRESSINGS




KERRACEL™ Ag Gelling Fiber Dressing

KERRACEL™ Ag Dressing is a gelling fiber silver dressing that uses the combination of an enhanced carboxymethyl cellulose (CMC) dressing and the advanced silver technology – Ag OXYSALTS™ Technology – to create the optimal wound healing environment.

KERRACEL™ Ag Gelling Fiber Dressing

Item code	Size	Eaches/Box	HCPCS code
CWL1157	2" x 2"	10 ea/bx	A6196
CWL1158	4" x 5"	10 ea/bx	A6197
CWL1159	6" x 6"	5 ea/bx	A6197
CWL1160	8" x 12"	5 ea/bx	A6198
CWL1168	0.75" x 12" ribbon	5 ea/bx	A6199
CWL1162	0.75" x 18" ribbon	5 ea/bx	A6199
CWL1163	16" x 16"	5 ea/bx	---




SILVERCEL™ NON-ADHERENT Antimicrobial Alginate Dressing

SILVERCEL™ NON-ADHERENT Antimicrobial Alginate Dressing with EASYLIFT™ Precision Film Technology is a non-woven pad composed of alginate, carboxymethylcellulose (CMC) and silver coated nylon fibers, with a non-adherent wound contact layer.

SILVERCEL™ NON-ADHERENT Antimicrobial Alginate Dressing

Item code	Size	Eaches/Carton/Box	HCPCS code
900202	2" x 2"	10 ea/ct - 5 ct/bx	A6196
900404	4 1/4" x 4 1/4"	10 ea/ct - 5 ct/bx	A6197
900408	4" x 8"	5 ea/ct - 5 ct/bx	A6197
900112	1" x 12"	5 ea/ct - 5 ct/bx	A6199



KERRACONTACT™ Ag Dressing

KERRACONTACT™ Ag Dressing is a silver dressing for chronic and acute wounds that utilizes Ag OXYSALTS™ Technology. Ag OXYSALTS™ Technology is a unique type of silver that is up to 6 times more powerful than Ag¹⁺.³

KERRACONTACT™ Ag Dressing

Item code	Size	Eaches/Box	HCPCS code
PN-09-0102	2" x 2"	10 ea/bx	A6206
PN-09-0103	4" x 5"	10 ea/bx	A6207
PN-09-0104	4" x 8"	10 ea/bx	A6207
PN-09-0105	6" x 6"	10 ea/bx	A6207
PN-09-0106	6" x 9"	10 ea/bx	A6208
PN-09-0108	8" x 16"	10 ea/bx	A6208
PN-09-0116	16" x 16"	5 ea/bx	A6208
PN-09-0120	4" x 48"	5 ea/bx	A6208
PN-09-0132	16" x 32"	5 ea/bx	A6208

PROMOTE GRANULATION




PROMOGRAN™ Matrix Wound Dressing

To help promote a moist wound environment conducive to granulation tissue formation, epithelialization, and rapid wound healing. PROMOGRAN™ Matrix is composed of collagen and Oxidized Regenerated Cellulose (ORC).



PROMOGRAN PRISMA™ Matrix

PROMOGRAN PRISMA™ Matrix contains collagen, ORC, and silver-ORC and has nominally twice the amount of Collagen/ORC material as the same size PROMOGRAN™ Matrix. Forms a soft, conformable, biodegradable gel.



FIBRACOL™ Plus Collagen Wound Dressing with Alginate

A combination of natural biopolymers combines the structural support of collagen and the gel forming properties of alginates into a sterile, soft, absorbent, and conformable topical wound dressing.



PROMOGRAN™ Matrix

Item code	Size	Eaches Carton/Box	HCPCS code
PG004	4.34 sq. in. hexagon	10 ea/ct - 4 ct/bx	A6021
PG019	19.07 sq. in. hexagon	10 ea/ct - 4 ct/bx	A6022

PROMOGRAN PRISMA™ Matrix

Item code	Size	Eaches Carton/Box	HCPCS code
MA028	4.34 sq. in. hexagon	10 ea/ct - 4 ct/bx	A6021
MA123	19.07 sq. in. hexagon	10 ea/ct - 4 ct/bx	A6022
MA032	3/8" x 3/8" x 12 5/8" rope	6 ea/ct - 4 ct/bx	A6024

FIBRACOL™ Plus Dressing

Item code	Size	Eaches Carton/Box	HCPCS code
2981	2" x 2"	12 ea/ct - 6 ct/bx	A6021
2982	4" x 4 3/8"	12 ea/ct - 6 ct/bx	A6022
2983	4" x 8 3/4"	6 ea/ct - 6 ct/bx	A6022
2984	3/8" x 3/8" x 15 3/4" rope	6 ea/ct - 6 ct/bx	A6024

EXUDATE MANAGEMENT



KERRAMAX CARE™ Dressing Family

A super-absorbent dressing that combines a unique wicking layer to distribute exudate throughout the dressing and absorptive polymers that absorb and retain high levels of exudate into a durable yet conformable dressing



KERRAFOAM™ Gentle Border Foam Dressing

Contains an absorbent core and foam pad to absorb and lock in exudate. The soft silicone adhesive contact layer helps to prevent painful dressing changes and skin stripping making this dressing ideal for both wound care and skin protection.



KERRAMAX CARE™ Dressing Family

	Item code	Size	Eaches/Box	HCPCS code
Non-Border	PRD500-025	2" x 2"	10 ea/bx	A6196
	PRD500-050	4" x 4"	10 ea/bx	A6196
	PRD500-100	5" x 6"	10 ea/bx	A6197
	PRD500-120	4" x 9"	10 ea/bx	A6197
	PRD500-240	8" x 9"	10 ea/bx	A6197
	PRD500-380	8" x 12"	5 ea/bx	A6198
	PRD500-600	8" x 20"	5 ea/bx	A6198
	PRD500-300	8" x 9" multisite	5 ea/bx	A6197
Silicone Border	PRD500-1174	4" x 4"	10 ea/bx	A6196
	PRD500-1175	6" x 6"	10 ea/bx	A6196
	PRD500-1176	8" x 8"	10 ea/bx	A6197
	PRD500-1177	6" x 10"	10 ea/bx	A6197

KERRAFOAM™ Gentle Border Foam Dressing



Item code	Size	Eaches/Box	HCPCS code
CWL1010	3" x 3"	10 ea/bx	A6212
CWL1011	4" x 4"	10 ea/bx	A6212
CWL1013	5" x 5"	10 ea/bx	A6212
CWL1132	6" x 6"	10 ea/bx	A6213
CWL1135	7" x 8" oval	10 ea/bx	A6213
CWL1012	4" x 8"	10 ea/bx	A6212
CWL1133	4" x 10"	10 ea/bx	A6213
CWL1134	4" x 12"	10 ea/bx	A6213
CWL1023	6.7" x 6.9" small sacrum	5 ea/bx	A6212
CWL1040	9" x 10" large sacrum	5 ea/bx	A6213
CWL1024	9.8" x 9" heel	5 ea/bx	A6213

EXUDATE MANAGEMENT




KERRAFOAM™ Simple Non Border Dressing

Designed to absorb excess exudate from the wound bed, while protecting surrounding skin to aid the healing of acute wounds. These non-adherent foam dressings are able to be cut to size.



TIELLE™ Non Adhesive Hydropolymer Dressing

TIELLE™ Non Adhesive Dressing with LIQUALOCK™ Technology is a cuttable hydropolymer foam dressing for low to heavily exuding wounds. The dressing contains LIQUALOCK™ Advanced absorption technology and uses a combination of MVTR and fluid absorption to manage exudate.




KERRALITE COOL™ Hydrogel Dressing

A soothing, debriding, and moisturizing dressing that provides the ideal environment for managing dry to lightly exuding sloughy wounds.




KERRACEL™ Gelling Fiber Dressing

KERRACEL™ is a 100% carboxymethyl cellulose primary dressing for use on a variety of exuding wounds. Soft and conformable, KERRACEL™ has a high absorbency and retains its integrity when removed from the wound or cavity.⁴



NU-DERM™ Hydrocolloid Wound Dressing Family

NU-DERM™ Hydrocolloid Dressings are designed to help maintain a moist wound environment. There are three variants in the range, available in a variety of shapes and sizes. All dressings are waterproof and will remain in place during showering.



NU-DERM™ Alginate Wound Dressing

NU-DERM™ Alginate Wound Dressing is a sterile, non-woven pad consisting of high G (guluronic acid) alginate and carboxymethylcellulose (CMC) fiber. NU-DERM™ Alginate Wound Dressing is designed to be fast gelling and to allow removal intact.



KERRAFOAM™ Simple Non Border Dressing

Item code	Size	Eaches/Box	HCPCS code
CWL1204	2" x 2"	10 ea/bx	A6209
CWL1206	4" x 4.3"	10 ea/bx	A6210
CWL1208	4" x 8"	10 ea/bx	A6210
CWL1207	6" x 6"	10 ea/bx	A6210
CWL1209	8" x 8"	10 ea/bx	A6211

TIELLE™ Non Adhesive Hydropolymer Dressing

Item code	Size (pad size)	Eaches carton/Box	HCPCS code
TLN0505EN	2" x 2"	10 ea/ct – 5 ct/bx	A6209
TLN1010EN	4" x 4"	10 ea/ct – 5 ct/bx	A6209
TLN1515EN	6" x 6"	10 ea/ct – 5 ct/bx	A6210
TLN1717EN	7" x 7"	5 ea/ct – 5 ct/bx	A6211
TLN2122EN	8 1/4" x 8 5/8"	5 ea/ct – 5 ct/bx	A6211

KERRALITE COOL™ Hydrogel Dressings

	KERRALITE COOL™ Hydrogel Dressings			
	Item code	Size	Eaches/Box	HCPCS code
Non-Adhesive	CWL1004	2.4" x 2.4"	5 ea/bx	A6234
	CWL1005	5" x 3.5"	5 ea/bx	A6235
	CWL1006	7" x 5"	5 ea/bx	A6235
Adhesive Border	CWL1007	3" x 3"	5 ea/bx	A6237
	CWL1008	4" x 4"	5 ea/bx	A6237
	CWL1009	6" x 6"	5 ea/bx	A6238

KERRACEL™ Gelling Fiber Dressing

Item code	Size	Eaches/Box	HCPCS code
CWL1032	2" x 2"	10 ea/bx	A6196
CWL1033	4"x 5"	10 ea/bx	A6197
CWL1034	6" x 6"	5 ea/bx	A6197
CWL1035	1" x 18" ribbon	5 ea/bx	A6199


NU-DERM™ Hydrocolloid Wound Dressing Family

Item code	Size	Eaches Tray/Case	HCPCS code
NU-DERM™ Thin Hydrocolloid Dressing			
HCT101	4" x 4"	10 ea/tray - 10 tray/cs	A6234
NU-DERM™ Standard Hydrocolloid Dressing			
HCF204	4" x 4"	5 ea/tray - 10 tray/cs	A6234
HCF208	8" x 8"	5 ea/tray - 4 tray/cs	A6236
NU-DERM™ Border Hydrocolloid Dressing			
HCB102	2" x 2"	20 ea/tray- 5 tray/cs	A6234
HCB204	4" x 4"	10 ea/tray- 16 tray/cs	A6234
HCB106	6" x 6"	5 ea/tray- 4 tray/cs	A6235
NU-DERM™ Sacral Hydrocolloid Dressing			
HCS100	6" X 7"	5 ea/tray - 4 tray/cs	A6237

NU-DERM™ Alginate Wound Dressing

Item code	Size	Eaches Carton/Box	HCPCS code
AWD112	1" x 12"	5 ea/ct - 5 ct/bx	A6199
AWD202	2" x 2"	10 ea/ct - 5 ct/bx	A6196
AWD404	4" x 4"	10 ea/ct - 5 ct/bx	A6196
AWD408	4" x 8"	5 ea/ct - 5 ct/bx	A6197

CONTACT LAYER DRESSING




ADAPTIC TOUCH™ Non-Adhering Silicone Dressing
ADAPTIC TOUCH™ Non-Adhering Silicone Dressing is a wound contact layer that contains soft tack silicone that assists with dressing application, conformability, and atraumatic removal.



ADAPTIC™ Digit Non-Adhering Dressing
ADAPTIC™ Digit Dressing is made of a unique, tacky, silicone-coated fabric covered with a tubular bandage. It is specifically designed for dressing fingers and toes.



ADAPTIC™ Non-Adhering Dressing
ADAPTIC™ Dressing is a primary dressing made of knitted cellulose acetate fabric and impregnated with a specially formulated petrolatum emulsion.



ADAPTIC TOUCH™ Dressing			
Item code	Size	Eaches Carton/Box	HCPSC code
500501	2" x 3"	10 ea/ct - 5 ct/bx	A6206
500502	3" x 4 1/4"	10 ea/ct - 5 ct/bx	A6206
500503	5" x 6"	10 ea/ct - 5 ct/bx	A6207
500504	8" x 12 3/4"	5 ct/bx	A6208


ADAPTIC™ Digit Non-Adhering Dressing		
Item code	Size	Eaches/Box
MAD003U	Small 3/4"	10 ea/bx
MAD013U	Med 1"	10 ea/bx
MAD023U	Large 1 1/8"	10 ea/bx
MAD042U	XL 1 1/4"	10 ea/bx
MAD062U	Large 1 1/8" toe	10 ea/bx

ADAPTIC™ Non-Adhering Dressing			
Item code	Size	Eaches Carton/Box	HCPSC code
2012	3" x 3"	50 ea/ct - 12 ct/bx	A6222
2013	3" x 8"	36 ea/ct - 6 ct/bx	A6223
2014	3" x 16"	36 ea/ct - 6 ct/bx	A6223
2015	3" x 8"	24 ea/ct - 6 ct/bx	A6223
2018	3" x 60" roll	10 rolls/bx	A6266
2019	5" x 9" sheet	12 ea/ct - 6 ct/bx	A6223

TRANSPARENT FILM




BIOCLUSIVE™ Plus Transparent Film Dressing
BIOCLUSIVE™ Plus Dressing is a self-adhesive dressing consisting of a transparent semi-permeable polyurethane film coated with an acrylic adhesive.




BIOCLUSIVE™ Plus Transparent Film Dressing			
Item code	Size	Eaches Carton/Box	HCPSC code
BIP0607	2 3/8" x 2 3/4"	10 ea/ct - 10 ct/bx - 15 bx/cs	A6257
BIP1012	4" x 4 3/4"	10 ea/ct - 10 ct/bx - 6 bx/cs	A6258
BIP1520	5 7/8" x 7 7/8"	10 ea/ct - 5 ct/bx - 6 bx/cs	A6258
BIP2030	7 7/8" x 11 3/4"	10 ea/ct - 5 ct/bx - 4 bx/cs	A6259

ODOR MANAGEMENT



ACTISORB™ Silver 220 Antimicrobial Binding Dressing
ACTISORB™ Silver 220 Antimicrobial Binding Dressing is an activated charcoal dressing with silver, enclosed in a non-adherent nylon sleeve.



ACTISORB™ Silver 220 Antimicrobial Binding Dressing			
Item code	Size	Eaches carton/box	HCPSC code
650220	2.6" x 3.8"	10 ea/ct - 5 ea/bx	A4649
105220	4.2" x 4.2"	10 ea/ct - 5 ea/bx	A4649
190220	4.2" x 7.6"	10 ea/ct - 5 ea/bx	A4649

COMPRESSION SYSTEM



DYNA-FLEX™ Multi-Layer Compression System
DYNA-FLEX™ Compression System is a multi-layer compression system that provides sustained compression for the management of venous leg ulcers and related conditions.

DYNA-FLEX™ Multi-Layer Compression System		
Item code	Eaches/Case	HCPSC Code
7025	8 units/cs	A6223- ADAPTIC™ Non-Adhering Dressing
		A6421 - 1st Layer
		A6430 - 2nd Layer
		A643B - 3rd Layer

For additional product details, please call **800-275-4524** or visit **myKCI.com**

References:

1. Sirelkhatim A, et. al. Review on Zinc Oxide Nanoparticles: Antibacterial Activity & Toxicity Mechanism. Nano-Micro Lett. 2015. 7(3):219-242. 2. Nolan, K. and Murmur, E. Moisturizers: Reality and the Skin benefits. Dermatologic Therapy. 2012;(25):229-233. 3. Spina C. Silver I, II, III: Chemical characteristics, properties, and anti-microbial activity. Crawford Healthcare Ltd. 4. Crawford Healthcare. CHC-R668. In vitro tests to analyze KERRACEL™ Dressing

Copyright 2019 KCI Licensing, Inc. All rights reserved. Unless otherwise noted, all trademarks designated herein are owned or licensed by KCI Licensing, Inc., KCI USA, Inc., Systagenix Wound Management, Ltd., or Crawford Healthcare, Ltd. LIT#34-L-033 • PRA-PM-US-01119 (01/19)

